

CALF IMPLANTS




SEBBIN PARIS
ESTHETIQUE & RECONSTRUCTION

Recreating harmonious bodies for a new vision of rebirth.

Reference	Length (mm)	Width (mm)	Projection (mm)	Volume (mL)
LS03120	120	40	15	40
LS03175	175	35	20	60
LS03180	180	42	21	85
LS03185	185	46	23	110
LS03190	190	50	25	140
LS03195	195	57	27	165

The dimensions of our implants can differ slightly from the data below. Production upon request of a textured envelope (reference LST).



- Created for reconstructive and aesthetic surgery.
- Insert pocket allowing easier insertion.



Pocket easing insertion

Our calf implants are sold sterile, single-use, by unit.

CONTINUITY Production know-how mastered for almost 30 years.

RELIABILITY Visual controls on 100% of implants at each step of the manufacturing process. Biological evaluation of end products.

TRACKING Each implant is identified by a unique serial number allowing tracking from the raw materials used to its implantation.

PEACE OF MIND Highly qualified personnel working in clean room.

This document is dedicated to health care professionals. According to the Medical Device Directive 93/42/EEC, the calf implants are Class IIb medical devices manufactured by Groupe SEBBIN and are intended to be used in the context of a plastic, reconstructive and cosmetic surgery. These devices are CE marked by the notified body mdc medical device certification GmbH number 0483. Please read carefully the instructions for use before using.

MANUFACTURER:



GROUPE SEBBIN SAS
 39, Parc d'Activités des Quatre Vents
 95650 Boissy l'Aillierie - France
 Tel.: +33 1 34 42 13 28
 Fax: +33 1 34 42 16 88
 www.sebbin.com

INTERNATIONAL OPERATIONS:

Sebbin Benelux SPRL:
 Tel.: +32 3 451 26 81 - Fax: +32 3 451 26 89
 Email: info@laboratoires-sebbin.be

Sebbin Deutschland GmbH:
 Tel.: +49 2102 70 88 00 0 - Fax: +49 2102 70 88 00 99
 Email: service@sebbin-implantate.de

Groupe Sebbin SAS, France:
 Tel.: +33 1 34 42 13 28 - Fax: +33 1 34 42 16 88
 Email: serviceclients@laboratoires-sebbin.fr

Sebbin Iberica Productos Medicos, S.L.:
 Tel.: +34 93 159 3003 - Fax: +34 93 676 3053
 Email: info@sebbin.es

Sebbin Suisse Sarl:
 Tel.: +41 43 336 72 30 - Fax: +41 43 336 72 36
 Email: info@sebbin.ch

Sebbin UK Ltd:
 Tel.: 0800 652 4028 - Fax: +4433 0088 9575
 Email: custserv@sebbin.uk