

FACIAL IMPLANTS




SEBBIN PARIS
ESTHETIQUE & RECONSTRUCTION

Recreating harmonious bodies for a new vision of rebirth.

Indications



Facial implants are indicated for aesthetic and reconstructive or corrective surgery. Different implants are offered:

- Short chin implants are designed for aesthetic improvement of the chin.
 - Long chin implants are designed for the same correction as standard chin implants, the difference consisting in an upper expansion of the lateral parts of the chin.
 - Round or anatomical malar implants are designed for aesthetic improvement of the cheekbones.
- 

Composition




Our facial implants are manufactured with a biocompatible silicone elastomer monobloc of medical quality, with long-term implantability. Their perforations allow the implant to be fixed to the tissues to favour its proper implementation.

Advantages




The material used is sufficiently supple and flexible to facilitate insertion via a small incision; it is sufficiently rigid to simulate the part of the face that is missing, or in need of correction. Facial implants keep their original shape and are not resorbed.




Short chin implant

Reference	Length (mm)	Width (mm)	Projection (mm)
LS 02 101	34	15	6
LS 02 102	44	16	7
LS 02 103	54	17	8




Long chin implant

Reference	Length (mm)	Small length (mm)	Width (mm)	Projection (mm)
LS 02 111	63	39	15	8
LS 02 112	73	49	15	9




Round malar implant

Reference	Length (mm)	Width (mm)	Projection (mm)
LS 02 221	38	20	4
LS 02 222	39	21	5
LS 02 223	41	23	5



Anatomical malar implant

Reference		Length (mm)	Width (mm)	Projection (mm)
Left	Right			
LS 02G 231	LS 02D 231	48	27	3
LS 02G 232	LS 02D 232	48	27	4
LS 02G 233	LS 02D 233	48	27	6
LS 02G 234	LS 02D 234	52	29	3
LS 02G 235	LS 02D 235	52	29	4
LS 02G 236	LS 02D 236	52	29	6



Our facial implants are sold sterile, single-use, by unit.

Our commitments

CONTINUITY

Production know-how mastered for almost 30 years.

RELIABILITY

Visual controls on 100% of implants at each step of the manufacturing process.
Biological evaluation of end products.

TRACKING

Each implant is identified by a unique serial number allowing tracking from the raw materials used to its implantation.

PEACE OF MIND

Highly qualified personnel working in clean room.

This document is dedicated to health care professionals. According to the Medical Device Directive 93/42/EEC, the facial implants are Class IIb medical devices manufactured by Groupe SEBBIN and are intended to be used in the context of a plastic, reconstructive and cosmetic surgery. These devices are CE marked by the notified body mdc medical device certification GmbH number 0483. Please read carefully the instructions for use before using.

MANUFACTURER:



GROUPE SEBBIN SAS
39, Parc d'Activités des Quatre Vents
95650 Boissy l'Aillierie - France
Tel.: +33 1 34 42 13 28
Fax: +33 1 34 42 16 88
www.sebbin.com

INTERNATIONAL OPERATIONS:

Sebbin Benelux SPRL:
Tel.: +32 3 451 26 81 - Fax: +32 3 451 26 89
Email: info@laboratoires-sebbin.be
Sebbin Deutschland GmbH:
Tel.: +49 2102 70 88 00 0 - Fax: +49 2102 70 88 00 99
Email: service@sebbin-implantate.de
Groupe Sebbin SAS, France:
Tel.: +33 1 34 42 13 28 - Fax: +33 1 34 42 16 88
Email: serviceclients@laboratoires-sebbin.fr
Sebbin Iberica Productos Medicos, S.L.:
Tel.: +34 93 159 3003 - Fax: +34 93 676 3053
Email: info@sebbin.es
Sebbin Suisse Sarl:
Tel.: +41 43 336 72 30 - Fax: +41 43 336 72 36
Email: info@sebbin.ch
Sebbin UK Ltd:
Tel.: 0800 652 4028 - Fax: +4433 0088 9575
Email: custserv@sebbin.uk