


Implant Mammaire Rond
Round Breast Implant



IMPLANTS CLASSIQUES | Naturgel™ Cohésif CLASSIC IMPLANTS | Cohesive Naturgel™

IMPLANTS FERMES | Naturgel™ Haute Cohésion FIRM IMPLANTS | High Cohesive Naturgel™

TEXTURES
TEXTURED

MICRO-TEXTURES
MICRO-TEXTURED

LISSES
SMOOTH

Profil Semi-moderé Semi-moderate profile	Profil Moderé Moderate profile	Profil Haut High profile																																																																																																																																																																										
N.A.	<table border="1"> <thead> <tr> <th colspan="5">LS 90</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 90 220</td><td>108</td><td>35</td><td>75</td><td></td></tr> <tr><td>LS 90 240</td><td>110</td><td>36</td><td>77</td><td></td></tr> <tr><td>LS 90 260</td><td>112</td><td>38</td><td>78</td><td></td></tr> <tr><td>LS 90 280</td><td>116</td><td>39</td><td>80</td><td></td></tr> <tr><td>LS 90 300</td><td>119</td><td>40</td><td>82</td><td></td></tr> <tr><td>LS 90 330</td><td>123</td><td>41</td><td>84</td><td></td></tr> <tr><td>LS 90 360</td><td>126</td><td>42</td><td>87</td><td></td></tr> <tr><td>LS 90 390</td><td>128</td><td>44</td><td>89</td><td></td></tr> <tr><td>LS 90 430</td><td>133</td><td>45</td><td>93</td><td></td></tr> <tr><td>LS 90 470</td><td>137</td><td>47</td><td>96</td><td></td></tr> </table>	LS 90					Ref.	Base	Proj.	"LVC"		LS 90 220	108	35	75		LS 90 240	110	36	77		LS 90 260	112	38	78		LS 90 280	116	39	80		LS 90 300	119	40	82		LS 90 330	123	41	84		LS 90 360	126	42	87		LS 90 390	128	44	89		LS 90 430	133	45	93		LS 90 470	137	47	96		<table border="1"> <thead> <tr> <th colspan="5">LS 91</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 91 225</td><td>97</td><td>46</td><td>73</td><td></td></tr> <tr><td>LS 91 245</td><td>99</td><td>47</td><td>75</td><td></td></tr> <tr><td>LS 91 265</td><td>101</td><td>49</td><td>78</td><td></td></tr> <tr><td>LS 91 285</td><td>103</td><td>50</td><td>80</td><td></td></tr> <tr><td>LS 91 305</td><td>105</td><td>52</td><td>82</td><td></td></tr> <tr><td>LS 91 335</td><td>109</td><td>54</td><td>85</td><td></td></tr> <tr><td>LS 91 365</td><td>112</td><td>56</td><td>87</td><td></td></tr> <tr><td>LS 91 395</td><td>114</td><td>57</td><td>89</td><td></td></tr> <tr><td>LS 91 435</td><td>117</td><td>60</td><td>91</td><td></td></tr> <tr><td>LS 91 475</td><td>122</td><td>61</td><td>95</td><td></td></tr> </table>	LS 91					Ref.	Base	Proj.	"LVC"		LS 91 225	97	46	73		LS 91 245	99	47	75		LS 91 265	101	49	78		LS 91 285	103	50	80		LS 91 305	105	52	82		LS 91 335	109	54	85		LS 91 365	112	56	87		LS 91 395	114	57	89		LS 91 435	117	60	91		LS 91 475	122	61	95																																																			
LS 90																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 90 220	108	35	75																																																																																																																																																																									
LS 90 240	110	36	77																																																																																																																																																																									
LS 90 260	112	38	78																																																																																																																																																																									
LS 90 280	116	39	80																																																																																																																																																																									
LS 90 300	119	40	82																																																																																																																																																																									
LS 90 330	123	41	84																																																																																																																																																																									
LS 90 360	126	42	87																																																																																																																																																																									
LS 90 390	128	44	89																																																																																																																																																																									
LS 90 430	133	45	93																																																																																																																																																																									
LS 90 470	137	47	96																																																																																																																																																																									
LS 91																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 91 225	97	46	73																																																																																																																																																																									
LS 91 245	99	47	75																																																																																																																																																																									
LS 91 265	101	49	78																																																																																																																																																																									
LS 91 285	103	50	80																																																																																																																																																																									
LS 91 305	105	52	82																																																																																																																																																																									
LS 91 335	109	54	85																																																																																																																																																																									
LS 91 365	112	56	87																																																																																																																																																																									
LS 91 395	114	57	89																																																																																																																																																																									
LS 91 435	117	60	91																																																																																																																																																																									
LS 91 475	122	61	95																																																																																																																																																																									
<table border="1"> <thead> <tr> <th colspan="5">LS 74</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 74 220</td><td>114</td><td>30</td><td>76</td><td></td></tr> <tr><td>LS 74 245</td><td>118</td><td>31</td><td>78</td><td></td></tr> <tr><td>LS 74 270</td><td>122</td><td>32</td><td>80</td><td></td></tr> <tr><td>LS 74 295</td><td>126</td><td>33</td><td>83</td><td></td></tr> <tr><td>LS 74 325</td><td>130</td><td>34</td><td>85</td><td></td></tr> <tr><td>LS 74 370</td><td>136</td><td>35</td><td>88</td><td></td></tr> <tr><td>LS 74 415</td><td>142</td><td>36</td><td>91</td><td></td></tr> <tr><td>LS 74 470</td><td>148</td><td>38</td><td>94</td><td></td></tr> </table>	LS 74					Ref.	Base	Proj.	"LVC"		LS 74 220	114	30	76		LS 74 245	118	31	78		LS 74 270	122	32	80		LS 74 295	126	33	83		LS 74 325	130	34	85		LS 74 370	136	35	88		LS 74 415	142	36	91		LS 74 470	148	38	94		<table border="1"> <thead> <tr> <th colspan="5">LS 70</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 70 220</td><td>108</td><td>35</td><td>75</td><td></td></tr> <tr><td>LS 70 240</td><td>110</td><td>36</td><td>77</td><td></td></tr> <tr><td>LS 70 260</td><td>112</td><td>38</td><td>78</td><td></td></tr> <tr><td>LS 70 280</td><td>116</td><td>39</td><td>80</td><td></td></tr> <tr><td>LS 70 300</td><td>119</td><td>40</td><td>82</td><td></td></tr> <tr><td>LS 70 330</td><td>123</td><td>41</td><td>84</td><td></td></tr> <tr><td>LS 70 360</td><td>126</td><td>42</td><td>87</td><td></td></tr> <tr><td>LS 70 390</td><td>128</td><td>44</td><td>89</td><td></td></tr> <tr><td>LS 70 430</td><td>133</td><td>45</td><td>93</td><td></td></tr> <tr><td>LS 70 470</td><td>137</td><td>47</td><td>96</td><td></td></tr> </table>	LS 70					Ref.	Base	Proj.	"LVC"		LS 70 220	108	35	75		LS 70 240	110	36	77		LS 70 260	112	38	78		LS 70 280	116	39	80		LS 70 300	119	40	82		LS 70 330	123	41	84		LS 70 360	126	42	87		LS 70 390	128	44	89		LS 70 430	133	45	93		LS 70 470	137	47	96		<table border="1"> <thead> <tr> <th colspan="5">LS 71</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 71 225</td><td>97</td><td>46</td><td>73</td><td></td></tr> <tr><td>LS 71 245</td><td>99</td><td>47</td><td>75</td><td></td></tr> <tr><td>LS 71 265</td><td>101</td><td>49</td><td>78</td><td></td></tr> <tr><td>LS 71 285</td><td>103</td><td>50</td><td>80</td><td></td></tr> <tr><td>LS 71 305</td><td>105</td><td>52</td><td>82</td><td></td></tr> <tr><td>LS 71 335</td><td>109</td><td>54</td><td>85</td><td></td></tr> <tr><td>LS 71 365</td><td>112</td><td>56</td><td>87</td><td></td></tr> <tr><td>LS 71 395</td><td>114</td><td>57</td><td>89</td><td></td></tr> <tr><td>LS 71 435</td><td>117</td><td>60</td><td>91</td><td></td></tr> <tr><td>LS 71 475</td><td>122</td><td>61</td><td>95</td><td></td></tr> </table>	LS 71					Ref.	Base	Proj.	"LVC"		LS 71 225	97	46	73		LS 71 245	99	47	75		LS 71 265	101	49	78		LS 71 285	103	50	80		LS 71 305	105	52	82		LS 71 335	109	54	85		LS 71 365	112	56	87		LS 71 395	114	57	89		LS 71 435	117	60	91		LS 71 475	122	61	95	
LS 74																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 74 220	114	30	76																																																																																																																																																																									
LS 74 245	118	31	78																																																																																																																																																																									
LS 74 270	122	32	80																																																																																																																																																																									
LS 74 295	126	33	83																																																																																																																																																																									
LS 74 325	130	34	85																																																																																																																																																																									
LS 74 370	136	35	88																																																																																																																																																																									
LS 74 415	142	36	91																																																																																																																																																																									
LS 74 470	148	38	94																																																																																																																																																																									
LS 70																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 70 220	108	35	75																																																																																																																																																																									
LS 70 240	110	36	77																																																																																																																																																																									
LS 70 260	112	38	78																																																																																																																																																																									
LS 70 280	116	39	80																																																																																																																																																																									
LS 70 300	119	40	82																																																																																																																																																																									
LS 70 330	123	41	84																																																																																																																																																																									
LS 70 360	126	42	87																																																																																																																																																																									
LS 70 390	128	44	89																																																																																																																																																																									
LS 70 430	133	45	93																																																																																																																																																																									
LS 70 470	137	47	96																																																																																																																																																																									
LS 71																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 71 225	97	46	73																																																																																																																																																																									
LS 71 245	99	47	75																																																																																																																																																																									
LS 71 265	101	49	78																																																																																																																																																																									
LS 71 285	103	50	80																																																																																																																																																																									
LS 71 305	105	52	82																																																																																																																																																																									
LS 71 335	109	54	85																																																																																																																																																																									
LS 71 365	112	56	87																																																																																																																																																																									
LS 71 395	114	57	89																																																																																																																																																																									
LS 71 435	117	60	91																																																																																																																																																																									
LS 71 475	122	61	95																																																																																																																																																																									
N.A.	<table border="1"> <thead> <tr> <th colspan="5">LS 50</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 50 220</td><td>108</td><td>35</td><td>75</td><td></td></tr> <tr><td>LS 50 240</td><td>110</td><td>36</td><td>77</td><td></td></tr> <tr><td>LS 50 260</td><td>112</td><td>38</td><td>78</td><td></td></tr> <tr><td>LS 50 280</td><td>116</td><td>39</td><td>80</td><td></td></tr> <tr><td>LS 50 300</td><td>119</td><td>40</td><td>82</td><td></td></tr> <tr><td>LS 50 330</td><td>123</td><td>41</td><td>84</td><td></td></tr> <tr><td>LS 50 360</td><td>126</td><td>42</td><td>87</td><td></td></tr> <tr><td>LS 50 390</td><td>128</td><td>44</td><td>89</td><td></td></tr> <tr><td>LS 50 430</td><td>133</td><td>45</td><td>93</td><td></td></tr> <tr><td>LS 50 470</td><td>137</td><td>47</td><td>96</td><td></td></tr> </table>	LS 50					Ref.	Base	Proj.	"LVC"		LS 50 220	108	35	75		LS 50 240	110	36	77		LS 50 260	112	38	78		LS 50 280	116	39	80		LS 50 300	119	40	82		LS 50 330	123	41	84		LS 50 360	126	42	87		LS 50 390	128	44	89		LS 50 430	133	45	93		LS 50 470	137	47	96		<table border="1"> <thead> <tr> <th colspan="5">LS 51</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LS 51 225</td><td>97</td><td>46</td><td>73</td><td></td></tr> <tr><td>LS 51 245</td><td>99</td><td>47</td><td>75</td><td></td></tr> <tr><td>LS 51 265</td><td>101</td><td>49</td><td>78</td><td></td></tr> <tr><td>LS 51 285</td><td>103</td><td>50</td><td>80</td><td></td></tr> <tr><td>LS 51 305</td><td>105</td><td>52</td><td>82</td><td></td></tr> <tr><td>LS 51 335</td><td>109</td><td>54</td><td>85</td><td></td></tr> <tr><td>LS 51 365</td><td>112</td><td>56</td><td>87</td><td></td></tr> <tr><td>LS 51 395</td><td>114</td><td>57</td><td>89</td><td></td></tr> <tr><td>LS 51 435</td><td>117</td><td>60</td><td>91</td><td></td></tr> <tr><td>LS 51 475</td><td>122</td><td>61</td><td>95</td><td></td></tr> </table>	LS 51					Ref.	Base	Proj.	"LVC"		LS 51 225	97	46	73		LS 51 245	99	47	75		LS 51 265	101	49	78		LS 51 285	103	50	80		LS 51 305	105	52	82		LS 51 335	109	54	85		LS 51 365	112	56	87		LS 51 395	114	57	89		LS 51 435	117	60	91		LS 51 475	122	61	95																																																			
LS 50																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 50 220	108	35	75																																																																																																																																																																									
LS 50 240	110	36	77																																																																																																																																																																									
LS 50 260	112	38	78																																																																																																																																																																									
LS 50 280	116	39	80																																																																																																																																																																									
LS 50 300	119	40	82																																																																																																																																																																									
LS 50 330	123	41	84																																																																																																																																																																									
LS 50 360	126	42	87																																																																																																																																																																									
LS 50 390	128	44	89																																																																																																																																																																									
LS 50 430	133	45	93																																																																																																																																																																									
LS 50 470	137	47	96																																																																																																																																																																									
LS 51																																																																																																																																																																												
Ref.	Base	Proj.	"LVC"																																																																																																																																																																									
LS 51 225	97	46	73																																																																																																																																																																									
LS 51 245	99	47	75																																																																																																																																																																									
LS 51 265	101	49	78																																																																																																																																																																									
LS 51 285	103	50	80																																																																																																																																																																									
LS 51 305	105	52	82																																																																																																																																																																									
LS 51 335	109	54	85																																																																																																																																																																									
LS 51 365	112	56	87																																																																																																																																																																									
LS 51 395	114	57	89																																																																																																																																																																									
LS 51 435	117	60	91																																																																																																																																																																									
LS 51 475	122	61	95																																																																																																																																																																									

Profil Semi-moderé Semi-moderate profile	Profil Moderé Moderate profile	Profil Haut High profile																																																																																																																																																																					
N.A.	<table border="1"> <thead> <tr> <th colspan="5">LSC 92</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 92 220</td><td>101</td><td>39</td><td>74</td><td></td></tr> <tr><td>LSC 92 240</td><td>105</td><td>39</td><td>76</td><td></td></tr> <tr><td>LSC 92 265</td><td>107</td><td>42</td><td>79</td><td></td></tr> <tr><td>LSC 92 285</td><td>109</td><td>43</td><td>80</td><td></td></tr> <tr><td>LSC 92 310</td><td>113</td><td>44</td><td>82</td><td></td></tr> <tr><td>LSC 92 330</td><td>117</td><td>45</td><td>84</td><td></td></tr> <tr><td>LSC 92 365</td><td>121</td><td>46</td><td>87</td><td></td></tr> <tr><td>LSC 92 395</td><td>124</td><td>48</td><td>89</td><td></td></tr> <tr><td>LSC 92 430</td><td>126</td><td>50</td><td>92</td><td></td></tr> <tr><td>LSC 92 475</td><td>130</td><td>51</td><td>95</td><td></td></tr> </table>	LSC 92					Ref.	Base	Proj.	"LVC"		LSC 92 220	101	39	74		LSC 92 240	105	39	76		LSC 92 265	107	42	79		LSC 92 285	109	43	80		LSC 92 310	113	44	82		LSC 92 330	117	45	84		LSC 92 365	121	46	87		LSC 92 395	124	48	89		LSC 92 430	126	50	92		LSC 92 475	130	51	95		<table border="1"> <thead> <tr> <th colspan="5">LSC 93</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 93 235</td><td>94</td><td>51</td><td>75</td><td></td></tr> <tr><td>LSC 93 255</td><td>96</td><td>53</td><td>78</td><td></td></tr> <tr><td>LSC 93 280</td><td>98</td><td>55</td><td>80</td><td></td></tr> <tr><td>LSC 93 300</td><td>101</td><td>56</td><td>82</td><td></td></tr> <tr><td>LSC 93 320</td><td>103</td><td>57</td><td>84</td><td></td></tr> <tr><td>LSC 93 350</td><td>106</td><td>58</td><td>86</td><td></td></tr> <tr><td>LSC 93 385</td><td>109</td><td>60</td><td>89</td><td></td></tr> <tr><td>LSC 93 415</td><td>112</td><td>62</td><td>91</td><td></td></tr> <tr><td>LSC 93 455</td><td>116</td><td>64</td><td>94</td><td></td></tr> </table>	LSC 93					Ref.	Base	Proj.	"LVC"		LSC 93 235	94	51	75		LSC 93 255	96	53	78		LSC 93 280	98	55	80		LSC 93 300	101	56	82		LSC 93 320	103	57	84		LSC 93 350	106	58	86		LSC 93 385	109	60	89		LSC 93 415	112	62	91		LSC 93 455	116	64	94																																																			
LSC 92																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 92 220	101	39	74																																																																																																																																																																				
LSC 92 240	105	39	76																																																																																																																																																																				
LSC 92 265	107	42	79																																																																																																																																																																				
LSC 92 285	109	43	80																																																																																																																																																																				
LSC 92 310	113	44	82																																																																																																																																																																				
LSC 92 330	117	45	84																																																																																																																																																																				
LSC 92 365	121	46	87																																																																																																																																																																				
LSC 92 395	124	48	89																																																																																																																																																																				
LSC 92 430	126	50	92																																																																																																																																																																				
LSC 92 475	130	51	95																																																																																																																																																																				
LSC 93																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 93 235	94	51	75																																																																																																																																																																				
LSC 93 255	96	53	78																																																																																																																																																																				
LSC 93 280	98	55	80																																																																																																																																																																				
LSC 93 300	101	56	82																																																																																																																																																																				
LSC 93 320	103	57	84																																																																																																																																																																				
LSC 93 350	106	58	86																																																																																																																																																																				
LSC 93 385	109	60	89																																																																																																																																																																				
LSC 93 415	112	62	91																																																																																																																																																																				
LSC 93 455	116	64	94																																																																																																																																																																				
<table border="1"> <thead> <tr> <th colspan="5">LSC 76</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 76 220</td><td>114</td><td>30</td><td>76</td><td></td></tr> <tr><td>LSC 76 245</td><td>118</td><td>31</td><td>78</td><td></td></tr> <tr><td>LSC 76 270</td><td>122</td><td>32</td><td>80</td><td></td></tr> <tr><td>LSC 76 295</td><td>126</td><td>33</td><td>83</td><td></td></tr> <tr><td>LSC 76 325</td><td>130</td><td>34</td><td>85</td><td></td></tr> <tr><td>LSC 76 370</td><td>136</td><td>35</td><td>88</td><td></td></tr> <tr><td>LSC 76 415</td><td>142</td><td>36</td><td>91</td><td></td></tr> <tr><td>LSC 76 470</td><td>148</td><td>38</td><td>94</td><td></td></tr> </table>	LSC 76					Ref.	Base	Proj.	"LVC"		LSC 76 220	114	30	76		LSC 76 245	118	31	78		LSC 76 270	122	32	80		LSC 76 295	126	33	83		LSC 76 325	130	34	85		LSC 76 370	136	35	88		LSC 76 415	142	36	91		LSC 76 470	148	38	94		<table border="1"> <thead> <tr> <th colspan="5">LSC 72</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 72 220</td><td>101</td><td>39</td><td>74</td><td></td></tr> <tr><td>LSC 72 240</td><td>105</td><td>39</td><td>76</td><td></td></tr> <tr><td>LSC 72 265</td><td>107</td><td>42</td><td>79</td><td></td></tr> <tr><td>LSC 72 285</td><td>109</td><td>43</td><td>80</td><td></td></tr> <tr><td>LSC 72 310</td><td>113</td><td>44</td><td>82</td><td></td></tr> <tr><td>LSC 72 330</td><td>117</td><td>45</td><td>84</td><td></td></tr> <tr><td>LSC 72 365</td><td>121</td><td>46</td><td>87</td><td></td></tr> <tr><td>LSC 72 395</td><td>124</td><td>48</td><td>89</td><td></td></tr> <tr><td>LSC 72 430</td><td>126</td><td>50</td><td>92</td><td></td></tr> <tr><td>LSC 72 475</td><td>130</td><td>51</td><td>95</td><td></td></tr> </table>	LSC 72					Ref.	Base	Proj.	"LVC"		LSC 72 220	101	39	74		LSC 72 240	105	39	76		LSC 72 265	107	42	79		LSC 72 285	109	43	80		LSC 72 310	113	44	82		LSC 72 330	117	45	84		LSC 72 365	121	46	87		LSC 72 395	124	48	89		LSC 72 430	126	50	92		LSC 72 475	130	51	95		<table border="1"> <thead> <tr> <th colspan="5">LSC 73</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 73 235</td><td>94</td><td>51</td><td>75</td><td></td></tr> <tr><td>LSC 73 255</td><td>96</td><td>53</td><td>78</td><td></td></tr> <tr><td>LSC 73 280</td><td>98</td><td>55</td><td>80</td><td></td></tr> <tr><td>LSC 73 300</td><td>101</td><td>56</td><td>82</td><td></td></tr> <tr><td>LSC 73 320</td><td>103</td><td>57</td><td>84</td><td></td></tr> <tr><td>LSC 73 350</td><td>106</td><td>58</td><td>86</td><td></td></tr> <tr><td>LSC 73 385</td><td>109</td><td>60</td><td>89</td><td></td></tr> <tr><td>LSC 73 415</td><td>112</td><td>62</td><td>91</td><td></td></tr> <tr><td>LSC 73 455</td><td>116</td><td>64</td><td>94</td><td></td></tr> </table>	LSC 73					Ref.	Base	Proj.	"LVC"		LSC 73 235	94	51	75		LSC 73 255	96	53	78		LSC 73 280	98	55	80		LSC 73 300	101	56	82		LSC 73 320	103	57	84		LSC 73 350	106	58	86		LSC 73 385	109	60	89		LSC 73 415	112	62	91		LSC 73 455	116	64	94	
LSC 76																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 76 220	114	30	76																																																																																																																																																																				
LSC 76 245	118	31	78																																																																																																																																																																				
LSC 76 270	122	32	80																																																																																																																																																																				
LSC 76 295	126	33	83																																																																																																																																																																				
LSC 76 325	130	34	85																																																																																																																																																																				
LSC 76 370	136	35	88																																																																																																																																																																				
LSC 76 415	142	36	91																																																																																																																																																																				
LSC 76 470	148	38	94																																																																																																																																																																				
LSC 72																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 72 220	101	39	74																																																																																																																																																																				
LSC 72 240	105	39	76																																																																																																																																																																				
LSC 72 265	107	42	79																																																																																																																																																																				
LSC 72 285	109	43	80																																																																																																																																																																				
LSC 72 310	113	44	82																																																																																																																																																																				
LSC 72 330	117	45	84																																																																																																																																																																				
LSC 72 365	121	46	87																																																																																																																																																																				
LSC 72 395	124	48	89																																																																																																																																																																				
LSC 72 430	126	50	92																																																																																																																																																																				
LSC 72 475	130	51	95																																																																																																																																																																				
LSC 73																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 73 235	94	51	75																																																																																																																																																																				
LSC 73 255	96	53	78																																																																																																																																																																				
LSC 73 280	98	55	80																																																																																																																																																																				
LSC 73 300	101	56	82																																																																																																																																																																				
LSC 73 320	103	57	84																																																																																																																																																																				
LSC 73 350	106	58	86																																																																																																																																																																				
LSC 73 385	109	60	89																																																																																																																																																																				
LSC 73 415	112	62	91																																																																																																																																																																				
LSC 73 455	116	64	94																																																																																																																																																																				
N.A.	<table border="1"> <thead> <tr> <th colspan="5">LSC 54</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 54 220</td><td>101</td><td>39</td><td>74</td><td></td></tr> <tr><td>LSC 54 240</td><td>105</td><td>39</td><td>76</td><td></td></tr> <tr><td>LSC 54 265</td><td>107</td><td>42</td><td>79</td><td></td></tr> <tr><td>LSC 54 285</td><td>109</td><td>43</td><td>80</td><td></td></tr> <tr><td>LSC 54 310</td><td>113</td><td>44</td><td>82</td><td></td></tr> <tr><td>LSC 54 330</td><td>117</td><td>45</td><td>84</td><td></td></tr> <tr><td>LSC 54 365</td><td>121</td><td>46</td><td>87</td><td></td></tr> <tr><td>LSC 54 395</td><td>124</td><td>48</td><td>89</td><td></td></tr> <tr><td>LSC 54 430</td><td>126</td><td>50</td><td>92</td><td></td></tr> <tr><td>LSC 54 475</td><td>130</td><td>51</td><td>95</td><td></td></tr> </table>	LSC 54					Ref.	Base	Proj.	"LVC"		LSC 54 220	101	39	74		LSC 54 240	105	39	76		LSC 54 265	107	42	79		LSC 54 285	109	43	80		LSC 54 310	113	44	82		LSC 54 330	117	45	84		LSC 54 365	121	46	87		LSC 54 395	124	48	89		LSC 54 430	126	50	92		LSC 54 475	130	51	95		<table border="1"> <thead> <tr> <th colspan="5">LSC 55</th> </tr> <tr> <th>Ref.</th> <th>Base</th> <th>Proj.</th> <th colspan="2">"LVC"</th> </tr> </thead> <tr><td>LSC 55 235</td><td>94</td><td>51</td><td>75</td><td></td></tr> <tr><td>LSC 55 255</td><td>96</td><td>53</td><td>78</td><td></td></tr> <tr><td>LSC 55 280</td><td>98</td><td>55</td><td>80</td><td></td></tr> <tr><td>LSC 55 300</td><td>101</td><td>56</td><td>82</td><td></td></tr> <tr><td>LSC 55 320</td><td>103</td><td>57</td><td>84</td><td></td></tr> <tr><td>LSC 55 350</td><td>106</td><td>58</td><td>86</td><td></td></tr> <tr><td>LSC 55 385</td><td>109</td><td>60</td><td>89</td><td></td></tr> <tr><td>LSC 55 415</td><td>112</td><td>62</td><td>91</td><td></td></tr> <tr><td>LSC 55 455</td><td>116</td><td>64</td><td>94</td><td></td></tr> </table>	LSC 55					Ref.	Base	Proj.	"LVC"		LSC 55 235	94	51	75		LSC 55 255	96	53	78		LSC 55 280	98	55	80		LSC 55 300	101	56	82		LSC 55 320	103	57	84		LSC 55 350	106	58	86		LSC 55 385	109	60	89		LSC 55 415	112	62	91		LSC 55 455	116	64	94																																																			
LSC 54																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 54 220	101	39	74																																																																																																																																																																				
LSC 54 240	105	39	76																																																																																																																																																																				
LSC 54 265	107	42	79																																																																																																																																																																				
LSC 54 285	109	43	80																																																																																																																																																																				
LSC 54 310	113	44	82																																																																																																																																																																				
LSC 54 330	117	45	84																																																																																																																																																																				
LSC 54 365	121	46	87																																																																																																																																																																				
LSC 54 395	124	48	89																																																																																																																																																																				
LSC 54 430	126	50	92																																																																																																																																																																				
LSC 54 475	130	51	95																																																																																																																																																																				
LSC 55																																																																																																																																																																							
Ref.	Base	Proj.	"LVC"																																																																																																																																																																				
LSC 55 235	94	51	75																																																																																																																																																																				
LSC 55 255	96	53	78																																																																																																																																																																				
LSC 55 280	98	55	80																																																																																																																																																																				
LSC 55 300	101	56	82																																																																																																																																																																				
LSC 55 320	103	57	84																																																																																																																																																																				
LSC 55 350	106	58	86																																																																																																																																																																				
LSC 55 385	109	60	89																																																																																																																																																																				
LSC 55 415	112	62	91																																																																																																																																																																				
LSC 55 455	116	64	94																																																																																																																																																																				

Dimensions en mm. Dimensions in mm.
Les dimensions indiquées dans ce présent document peuvent présenter une variation de 5%. The dimensions indicated into this document may vary of 5%.